

Formidabel

SIMPLY UNIQUE

FMB
care



Contents & Editorial

Actively and constantly improving the care situation for nursing home residents and hospital patients and, of course, for their carers – that is our aim at FMB care. And the most important building block in the achieving of this objective is our care bed series Formidabel. Day in, day out we work tirelessly in our Development department, in our production halls and with our partners, resolutely intent on the continued advancement of Formidabel. The result is that we have succeeded in presenting a care bed which is easy to operate and which offers a unique degree of comfort with a modern look that provides more cosiness. And into the bargain this bed is also amazingly cost-effective.

But even if our beds make us feel that we have achieved the ultimate in care sophistication, we do not rest on our laurels, on our solutions that most definitely set standards in our field of work. Fired up by ambition, we continue pursuing our company aim and strike out in new directions together with our customers. Not only together but as partners - that is the way we develop and produce beds, furniture and furnishing concepts with our clients, designs which can adapt to the requirements of every individual customer. Simply because we believe that a care bed should be as multifaceted and unique as Life itself. In other words, simply Formidabel - impressive and inspirational.

Contents & Editorial

Solutions for more freedom in the world of care

- A modular solution
- Side protection **Unique**
- Scenario handset
- Integrated bed extension
- Formidabel UltraLow
- Digital solutions

Formidabel – Care beds – the next level

- Formidabel **First**
- Formidabel **Pro**
- Formidabel **Eco**
- Formidabel **UltraLow**
- Side protection – **Unique, Solid Up, Solid, Basic**
- Side protection **Unique**
- Side protection **Solid Up**
- Side protection **Solid & Basic**

Equipment & Accessories

- Lying surfaces, Safer Exit and mattresses
- Castors, lifting poles, wall spacers and lamps
- Shapes and colours
- Bedside tables and overbed tables

care & interior

- The combination
- A pleasant ambience

Together we can achieve more

- Understanding
- Responsibility
- Legal notices / notes

„The market has changed.
It is now much more global, more challenging.
But this transformation does not
mean that our focus has changed.
For us the human aspect of our work still
steadfastly remains the top priority.“

MARTIN FERBER

Solutions for more freedom in the world of care

The person is always the most important factor in the equation – everyone who works in the care sector shares this conviction. So our task is quite unequivocal: the nurturing of innovative ideas and the creating of supreme products to make life easier and to improve its quality for everyone requiring care and for all those providing such care. Here at FMB we devote ourselves to this task wholeheartedly and strive towards its success every day.



A MODULAR SOLUTION






Absolutely unique variable

Every home resident, every hospital patient has special needs. Therefore we think it should always be possible to upgrade and convert a care bed along modular lines so that it can be flexibly and fuss-free adapted to suit individual requirements on the spot in each and every location. This was the thinking behind Formidabel and the reason why we developed it to become the most wide-ranging modular care bed on the market. The lying surface, the castors, the brake system or the side safety system or even the shape and colour of the headboard and the foot of the bed can all be combined in a constellation tailored to suit your ideas.

And this applies to all our beds, whether it's a First, Pro, Eco or an UltraLow. With FMB care you can rest assured that you will always get the care bed that you really need, the bed that will actively assist the caring work in your particular establishment. No compromises, totally versatile, absolutely innovative. In other words, simply Formidabel - impressive and inspirational.

We would love to talk to you personally and advise you in detail about all our modular solutions. Together we will be able to work out the special needs of your care home and then configure your bespoke Formidabel.

+ Advantages

-  Especially cost-effective
-  Unique and easy configuration
-  Powerful combinations
-  Solutions for more freedom
-  Always suitably equipped



SIDE PROTECTION UNIQUE

Freedom with safety

When developing our side protection system – and incidentally Unique was as clear-cut as it was ambitious – this was our aim: to reconcile the maximum amount of freedom with the maximum amount of safety. Spurred on by the belief that we could sustainably improve a care home resident's quality of life by means of a bed with innovative side protection, we went on to ultimately develop the perfect system. Fast, individual, safe. Unique.






Acting as a complement to the modular structure of our Formidabel series, the Unique side protection features adapt very easily to the needs of residents and to the requirements of hospital and nursing home operators. The maximum range of solutions, which can be used for speed, as a way of saving space and indeed in so many variable ways, constitutes a huge advantage for a carer team.

Unique enables the necessary protection for the resident while at the same time allowing the greatest amount of freedom possible. And this at all times. And thanks to their self-explanatory function, every Unique element can be mounted with very little effort and without the need for any tools. Scope and size are customised to suit protection requirements in the various situations.

Safety and transparency

Complete protection with four side protection elements, which can also be fixed just at the head end of the bed, the Unique side protection system with its transparent and very pleasant mesh fabric gives the user of the bed a particular feeling of freedom and also prevents any risk of injury caused by limbs being jammed. Thanks to its infinitely variable setting and its standard in accordance with the required stability factor, the side protection element also serves as a safe mobilisation aid for the user of the bed.

+ Advantages

-  Modular structure and consequently particularly cost-effective
-  No tools needed and easy to use
-  Innovative net structure provides more safety
-  Quick and easy cleaning
-  Maximum feeling of freedom for the user of the bed



Scenario 1



Scenario 2



Scenario 4



Scenario 3



SCENARIO HANDSET







Solutions with buttons

Good solutions operate smoothly. The Scenario handset makes the daily routine of residents and carers as easy as humanly possible: just by pressing a button.


Easy-to-understand symbols help the resident to intuitively have control over everyday situations just by pressing a button: scenarios such as finding a comfortable position to relax, a sleeping position or achieving the right height for the bed in order to get up easily and safely. The result is a very real potential lightening of the workload for carers. And for those needing care, it means more independence and consequently a higher level of social integration.

In addition to the programmable resident's mode, the configurable care mode and the lockdown mode work together to create a safe operating of the device. In short: we are talking more quality of life and more enjoyment at work.

+ Advantages

-  Programmable scenarios for more independence
-  Easy-to-understand symbols for intuitive operating
-  Button illumination provides excellent orientation even when it is dark
-  Functions can be selected and then blocked
-  Significant lightening of staff workload
-  Optional: Scenario handset for T16 / T17





BED EXTENSION WITH NO TOOLS REQUIRED

Longer and so much better




Every care home resident or hospital patient is different. So it is all the more important that a care bed can be flexibly adapted to suit the prevailing needs and conditions. And that has to happen with very little effort. The way in which the bed length of Formidabel First and UltraLow can be extended without the need for any tools is the perfect example of how good solutions can also be very simple.


Just a few quick moves is all that it takes

It is vital that a care bed can be adapted swiftly and in a straightforward manner to suit the height of each resident or individual patient. Thanks to the integrated bed extension feature, which comes

as a standard fitting, all our Formidabel models can be easily extended. With the First and UltraLow models that are available with the Unique side protection, the bed extension, carried out without the need for any tools, leads to even more flexibility and saves a great deal of time, too. By means of a practical release mechanism, these models can be quickly extended from 2 metres in length to 2.10 metres or even up to 2.20 metres. In this way it is possible to react at short notice if beds are to be used for different residents or patients with the result that the new “inhabitant” of the bed can be given the greatest comfort immediately.

+ Advantages

-  Extendible to 220 cm
– pure simplicity and no in-house technical work required
-  No tools needed, just a few moves of the hand do the job
-  Efficient procedures when bed becomes occupied by a new resident or patient



ULTRALOW




Safe proximity to the floor

Our UltraLow is the latest building block in the Formidabel series. This is where you will find all cutting-edge solutions and functions brought together in order to create the greatest sense of wellbeing possible for the care home resident. When the relevant handset scenario option is selected, the bed's lying surface is lowered to a height of just 15 cm away from the floor. Thus, UltraLow becomes instrumental in reducing consequences of any possible fall suffered by the resident. Furthermore, residents find it easier to get their bearings in the room, especially at night, because they can touch the floor with their hands and in this way feel much more confident of their position.

An ideal combination

The protection offered by the Unique side safety system is the ideal complement to the lowest lying position. On account of its infinitely variable adjustability, the side safety system serves as a safe aid for getting up out of bed or its handrail at mattress level can function as a cognitive signal warning of the risk of falling out of bed.

+ Advantages

-  Specific fall prevention by setting the lowest bed position
-  Unparalleled combination made up out of 15 cm proximity to the floor & the Unique side protection feature
-  Integrated fittings and planning reliability in every position







PRACTICAL SYNERGIES

Digital solutions

Two objectives are always of crucial importance for us: the improvement of care safety and quality, the lightening of the carers' workload. In combination with the Formidabel care bed, we can see how digital solutions can hugely support carers and enhance the quality of life for people in need of care. Our solution called Safer Exit, sensor-supported and developed together with the care sector, is an example of a technology that can be implemented simply and successfully. In nursing homes it is a fact that falls can quite often happen and the consequences thereof can be most serious. If a resident who is at particularly at risk from falling gets up from bed, the sensor system, suitably programmed for the individual resident, alerts the nursing staff via radio, thus enabling them to be on the spot quickly.

And especially when one considers challenges such as finding specialist staff and keeping them, our handset, with its selection of various scenarios typical of the care situation, offers a pragmatic, digital solution with invaluable benefits. Programmed functions, such as choosing the position for getting up out of bed and the comfortable position for relaxing or the operating of the underbed light, can all be selected just by the pressing of a button. These features intuitively make for an easier daily routine for the residents and they promote independence.

+ Advantages

-  More freedom and safety for high-risk residents prone to falling
-  Extra support at hand during night and weekend shifts
-  Fast reaction times in the event of falls
-  Retrofittable

“With a great amount of input and know-how, FMB puts into practice sustainable and innovative solutions as early as the conceptual stage. And because all the elements are optimally in tune with one another, the greatest potential is achieved.”

ANTHONY DONADINI

Moments Furniture, Switzerland

Formidabel – **Care beds** – the next level

The Formidabel care bed series is the result of more than 30 years of experience married together with the courage to be a trailblazer and with the understanding of problems and challenges found in the world of daily nursing care. It's time to find out what makes Formidabel simply unique and what it can uniquely do.



Formidabel First

The first of its kind

The care bed Formidabel First has made its mark on this sector, defining it anew by combining the highest possible degree of safety with the maximum amount of freedom. And it has done all this in an incredibly unique way. First and foremost because of its one-of-a-kind side protection system aptly named Unique, side safety that can be implemented along modular lines. Thanks to its continuous adjustment settings and state-of-the-art mesh fabric, a resident using this bed is restricted only as much as is necessary and therefore enjoys a feeling of freedom at all times.

Basic features

4 section lying surface
Metal slats, width 90 cm

Height adjustment
23 - 80 cm

Bed extension
Integrated, no tools needed

Double castors Ø 50 mm
With plastic cover

Handset
Scenario

Safe working load
225 kg / opt. 275 kg







Head position/ low position

2 + 2 brake system

Technology that makes Life easier

With the Unique side protection feature, you can react immediately to every care situation that arises. No special tools are required to fit this side safety system to the care bed. Its added bonus is that it also functions as a practical aid for the resident or patient to get up out of bed. Furthermore, because this bed can be extended also without the need for any tools, it is possible to adjust the Formidabel First to suit the requirements of every resident or patient using the bed. The handset allows the everyday scenarios to be selected as required and it ensures that the wide-ranging functions of this premium care bed, guaranteeing care and comfort for the user, are constantly available.

+ Compatible

-  **Side protection**
Unique (S. 30)
Solid Up (S. 32)
-  **Lying surface**
 - Special widths:
100 cm or 120 cm
 - Synthetic underlays,
Width 90 cm (p. 38)
-  **Covers**
Metal casing for
Ø 50 mm castors (p. 40)
-  **Double castors**
 - With Ø 75 mm and
2 + 2 brake system,
optional 4 castors
with a central brake
system (p. 40)
 - With Ø 100 mm and
2 + 2 brake system,
optional 4 castors
with a central brake
system (p. 40)
-  **Headboards and
footboards** (p. 43)
-  **Additional features**
(from page 38 onwards)



Formidabel Pro

Tried and tested taken to the next level

The Pro bed can be easily combined with the two side protection systems known as Solid and Basic. This bed symbolises traditional functions that chime together in harmony with brand new reliable solutions for the daily nursing routine.

Good sound care and workload reduction

In addition to a selection of various castor and brake systems, the care beds in the Pro series have on offer innumerable options when it comes to a choice of equipment. With the divided side protection known as Solid, which particularly facilitates getting in

and out of the bed thanks to its asymmetrical division of 40:60, this Formidabel model, too, provides safety without imposing measures that would deprive patients of their freedom. It also enables residents and patients to enjoy a lot of time devoted to valuable caring activities. The huge variety of functions such as the standard handset covering a selection of scenarios, the integrated bed extension facility, the high degree of stability, a safe working load and the kaleidoscope of design options perfectly round off the complete range of features and equipment offered by Formidabel Pro.

Basic features

4 section lying surface
Metal slats, width 90 cm

Height adjustment
25 - 80 cm

Bed extension
Integrated

Double castors Ø 50 mm
With plastic cover

Handset
Scenario

Safe working load
225 kg

Head position/ low position

2 + 2 brake system

+ Compatible

- + Side protection**
 - Solid (p. 35)
 - Basic (p. 34)
- + Lying surface**
 - Special width: 110 cm
 - Synthetic underlays, Width 90 cm (p. 38)
- + Covers**

Metal casing for
Ø 50 mm castors (p. 40)
- + Double castors**
 - With Ø 75 mm and 2 + 2 brake system, optional 4 castors with a central brake system (p. 40)
 - With Ø 100 mm and 2 + 2 brake system, optional 4 castors with a central brake system (p. 40)
- + Headboards and footboards** (p. 43)
- + Additional features**
(from page 38 onwards)



Formidabel

Eco

Kind to your budget

Formidabel Eco has been pared down to the most necessary and the most important functions, which are there to make life easier for residents and patients and for the nursing staff. The result is that the beds of the Eco series are particularly budget-friendly but nevertheless they do not have to go without the tried-and-tested stability and functions of a high-class care bed. As well as the side protection systems of Unique and Solid, the Eco model is also available with the more favourably-priced alternative of Basic.

Basic features

4 section lying surface
Metal slats, width 90 cm

Height adjustment
26 - 80 cm

Bed extension
Integrated

Double castors Ø 50 mm
With plastic cover

Handset
T7

Safe working load
175 kg

2 + 2 brake system

Individual consultation

We will be only too pleased to arrange a personal meeting so that we can give you detailed information on the qualities and benefits of the individual beds making up the Formidabel series. Together we will be able to find the Formidabel which is the perfect solution for your requirements.

+ Compatible

+ Side protection

- Unique (p. 30)
- Solid (p. 35)
- Basic (p. 34)

+ Lying surface

Synthetic underlays,
Width 90 cm (p. 38)

+ Double castors

- With Ø 100 mm and single wheel locking system (p. 40)
- With Ø 100 mm and 4 castors with a central brake system (p. 40)

+ Headboards and footboards (p. 43)

+ Additional features (from page 38 onwards)



Formidabel UltraLow

Fall prevention taken to the next level

With Formidabel UltraLow we bring together all our solutions for a cutting-edge care bed in a really unique way. A very low position close to the floor and our Unique side protection provide fantastic safety for residents and patients while still affording them freedom and independence.

The perfect addition

The lying surface of the UltraLow bed can be dropped down especially low to achieve a height of just 15 cm away from the floor. This means that with its safe proximity to the floor it contributes great-

ly to the prevention of falls without engaging any freedom-restraining measures. Together with the Unique side protection feature and its optimally coordinated modules, UltraLow boasts an extensive range of first-class fittings and equipment. Panels and side protection elements can be replaced without the need for any tools. On account of its infinitely variable adjustability, the Unique side protection feature serves as a safe aid for getting up out of bed or its handrail at mattress level can function as a cognitive signal warning of the risk of falling out of bed. This Formidabel model epitomises a high degree of both flexibility and also planning reliability for in-patient care.

Features

4 section lying surface
Metal slats, width 90 cm

Height adjustment
15 - 80 cm

Bed extension
Integrated, no tools needed

Double castors Ø 50 mm
With plastic cover

Handset
Scenario

Safe working load
225 kg

Head position/ low position

2 + 2 brake system

+ Optimally combined

- + Absolutely unique combination of optimally coordinated modules
- + For maximum freedom and the greatest possible safety
- + Wide selection of designs for a personalised and harmonious, cosy and welcoming atmosphere
- + **Side protection**
 - Unique (S. 30)
 - Solid Up (S. 32)
- + **Lying surface**
 - Special widths: 100 cm

“When it comes to health care, time is
an especially valuable treasure.
Our solutions guarantee that in the
hustle and bustle of daily care more time
remains for the most important thing,
the person being cared for.”




MARIUS FERBER

Side protection

Unique, Solid Up, Solid, Basic

It is no easy task to give a person in need of care the greatest possible amount of safety while maintaining the maximum amount of freedom for that person. But it is also a task which is monumentally important. And for that reason, on the one hand, we have optimised the old, tried-and-tested systems but, on the other, we have forged ahead to break new ground.

Side protection Unique



Panels

For models which need no side protection or safety rail system

Part protection

At the head end with two side protection elements for mobility and safety

Full protection

Provided by four side protection elements placed along the entire length of the bed



F P E UL

The Unique side protection system marries together maximum freedom, the greatest possible safety and luxurious comfort. The single elements are meticulously aligned with one another so as to enable countless combination options suitable for every situation and for every kind of requirement. Every Unique element can be mounted in exactly the place where it is really needed. And this can be done with just a few moves of the hand. No tools whatsoever are required. Offering complete protection provided by four side protection elements, which can also be fixed just at the head end of the bed, the Unique side protection system with its transparent,

pleasant and easy-to-clean mesh fabric gives users of the bed a particular feeling of freedom and wellbeing and also prevents any risk of injury caused by limbs being jammed. And thanks to its infinitely variable setting and its high level of stability, the protection element serves as a pleasant and safe aid for getting up out of bed and as a mobility aid. Drawn up to its maximum height, it can also function as a cognitive signal warning of the risk of falling out of bed or it can provide full protection against this. This variety of possibilities is what makes Unique live up to its name: it is so simply so unique in the world of care.

For Formidabel
First, UltraLow, Eco

Adjustability
Infinitely variable

Usable as follows
Modules freely selected as panels or as full protection

Multi-piece
Yes

Side protection
Unique

Compatibility

F Formidabel First

P Formidabel Pro


E Formidabel Eco

UL Formidabel UltraLow

Compatible


Not compatible

Side protection Solid Up




Panels

For models which need no side protection or safety rail system




Part protection

At the head end with two side protection elements for mobility and safety



Full protection

Provided by four side protection elements placed along the entire length of the bed



F P E UL

The side protection Solid Up ideally complements the product family of FMB care. Like Unique, Solid Up also offers - unlike conventional side guards - solutions that are tool-free and variable in just a few steps, if required. Depending on the care situation, the elements can be attached to the bed in different combinations. As side panels, if no side guard is needed, or as side guard elements for residents in need of special protection. A side guard attached on the head side supports entry and exit due to the asymmetrical division of the elements and thus offers additional safety. The maximum stability and flexibility of the Solid Up guarantee suitable, sustainable and budget-friendly solutions.

Adjustability

Two positions

Multi-piece

Yes

For Formidabel

First and UltraLow

Usable as follows

On both sides, on one side

Side protection Solid Up

Compatibility

F Formidabel **First**

P Formidabel **Pro**

E Formidabel **Eco**

UL Formidabel **UltraLow**

Compatible

Not compatible

Side protection Basic



F P E UL

The Basic side protection runs continuously around the side of the bed and, when fully pulled up with just a few moves of the hand, it provides protection and safety for those residents or patients who are particularly at risk of falling – after, of course, due consideration has been given to all legal conditions in line with FRM (freedom-restraining measures).

Adjustability
Three positions

For Formidabel
Pro and Eco

Multi-piece
No

Usable as follows
On both sides, on one side





In order to allow carers to react flexibly, swiftly and with the appropriate respect when protecting their residents and patients, our Solid side protection has a solution that meets such demands. It has been designed in such a way that its asymmetrical 40:60 division excellently facilitates getting into and out of the bed while at the same time offering the greatest possible safety as an aid for patients to rise. Because of all the various settings and possible adjustments, Solid guarantees stability and functionality when and where these are required.

Adjustability

Three positions

Multi-piece

Yes

For Formidabel

Pro and Eco

Usable as follows

On both sides, on one side

Side protection Solid

Compatibility

F Formidabel **First**

P Formidabel **Pro**

E Formidabel **Eco**

UL Formidabel **UltraLow**

Compatible

Not compatible

“We are constantly rethinking our care beds for the simple reason that different kinds of requirements demand personalised and optimally adjusted solutions.”

MARCEL KEUKER

FMB care, Netherlands

Equipment & **Accessories**

The right accessories mean that every kind of Formidabel bed can be easily and effortlessly tailored to suit the situation of the resident or patient concerned. Not only is that good for the comfort factor and the safety issues but it also means that use can be made of the available budgets in the best and most efficient way.



Lying surface with metal slats

The lying surface, divided into 4 sections with the retractable backrest for a pressure-relieving sitting position and with the integrated lying surface extension, creates a comfortable lying position second to none. All the lying surfaces of the Formidabel series are fitted with metal slats as a basic version. These slats are easy to clean and they help to create good ventilation of the bed.

Width
90 - 120 cm

Removable
No

Lying surface with synthetic underlays

Both the comfort and also the hygiene of the ergonomically divided up lying surface can be further optimised thanks to the easily removable synthetic underlays. The retractable backrest leading to relief of the pelvic area or the use of the integrated lying surface extension and also our handset offering the selection of various scenarios are all features that create supreme comfort and convenience for all living and working with such a bed.

Width
90 cm

Removable
Yes

Safer Exit (Illustration similar)

This sensor-supported solution can be fitted on to our lying surfaces quickly and in a totally straightforward way. Depending on the resident's mobility and risk of falling, the suitably programmed sensor system reports to the caring staff via radio that a patient has left the bed or possibly fallen out of bed. The carers can then react in an instant. Safer Exit does not only provide support for the daily caring routine but it also ramps up the level of safety enjoyed by residents and patients. And this is especially the case when it is used in conjunction with our beds that can be dropped to a really low position and with the Unique side protection system: it is truly an effective innovation that enables the reduction of freedom-restraining measures.

Mattresses

The variety of FMB care mattress models, ranging from standard through to comfort and also those for prevention of certain conditions, is especially tailored to suit the care environment. These mattresses offer residents and patients an exceptionally comfortable sleeping and lying position while at the same time fulfilling all the hygienic requirements demanded of the caring situation. Features such as foam composition, density, reinforced perimeter zones or a particular profile provide the appropriate structure, which supports not only mobility but which also brings about a pressure-relief effect.

For more details, please speak to our consultants. And we feel really sure that you will find the suitable mattress model for your establishment in our product portfolio.

Height	Width	Length
12 - 15 cm	90 - 120 cm	200 - 220 cm






Double castors Ø 50 mm

All double castors having a diameter of 50 mm are fitted with a matching plastic cover – and with First and Pro there is also the option of a metal cover.



Plastic cover

Removable cover for double castors having a diameter of 50 mm – for more homeliness and easy cleaning.



Metal cover

Removable cover for double castors having a diameter of 50 mm – both homely and hygienic.



Double castors Ø 75 mm

For easy and quiet moving and manoeuvrability of the bed, here with the practical 2 + 2 brake system.



Double castors Ø 100 mm

Super mobile exhibiting easy and quiet transport features, here with the convenient 4-castor central brake system.



Castors Ø 100 mm

Easy manoeuvrability and castor movement, here with single wheel locking system.



Wall spacer

Whether it is fixed at the head, foot and/or on the side of a Formidabel care bed, the result is always the same: you can rely on the wall spacer to protect the wall, the home's furnishings and the bed all in the same way.




Lifting pole with retracting strap

Height adjustment of the trapeze grab bar is achieved quickly and so flexibly by means of the retractable strap system that rolls up. In this way independence and mobility are both promoted.



Lifting pole with manual strap adjustment

The resident or patient can pull himself or herself up using the trapeze grab bar and then take up a new position more easily. The strap is manually adjusted at the suitable height.




Wall deflection rollers

Fitted to the headboards and footboards, these deflection rollers reliably protect furniture and furnishings against damage that could occur during the daily caring routine.



Culta bed reading lamp

The Culta LED reading lamp stands out because of the first-class materials that go into its making. The pliable swan neck and the rotating lamp head ensure the flexibility boasted by this lamp.



Handset holder

The handset is really easy to reach for residents, patients and carers thanks to the practical handset holder.

Compatibility

F Formidabel First

P Formidabel Pro

E Formidabel Eco

UL Formidabel UltraLow

Compatible

Not compatible

SHAPES AND COLOURS




As you like it

Good design is not just a question of looking good but it also enables aesthetics to chime together in harmony with functional properties. Modern wood décor bestows an attractive and high-quality look to our care furniture, creating an impression of luxury. Such decoration conveys a feeling of warmth and yet at the same time it complies with the high standards regarding cleaning and disinfection.

Your taste is the deciding factor

With the extensive selection of shapes and colours at our fingertips, together with you we can design every Formidabel care bed to be completely in tune with your taste and your décor. As well as our extensive range of ultra modern wood décor, we also have textile fabrics and upholstered headboards and footboards in low-maintenance fabrics – with no compromises regarding their use or function. All this provides countless creative options and ideas. The result is a cosy picture of homeliness felt by one and all.

+ Advantages

-  **Customised**
The huge selection of decorative possibilities on offer means that there are countless design options.
-  **Artistically expressive**
Textile fabrics, upholstered elements or special decorations all play their part in creating a distinctive kind of ambience.
-  **Practical**
All materials used can be quickly and easily cleaned and they all comply with hygienic requirements.

Compatible with side protection
Unique, Solid Up, Solid and Basic



F
P
E
UL

RD



F
P
E
UL

RG



F
P
E
UL

RO



**Décor
 Selection**
 Current décor and
 trends on offer
 can be found at
fmb-care.de

Compatible with side protection
Unique, Solid and Basic



F
P
E
UL

BG




F
P
E
UL

BO



F
P
E
UL

BN



Compatible with side protection
Unique




F
P
E
UL

RN



F
P
E
UL

PF



F
P
E
UL

PN



F
P
E
UL

BN /
 Covers

Compatibility

F Formidabel **First**

P Formidabel **Pro**

E Formidabel **Eco**

UL Formidabel **UltraLow**

Compatible

Not compatible

BEDSIDE TABLES AND OVERBED TABLES

The perfect companions

Independence is undeniably an important factor in the everyday life of residents and patients. A big part of this is being able to access personal storage space, shelves and surfaces. Our range of bedside tables and overbed tables is optimally tailored in design and function to suit our care beds and also to complement care needs.

The easy route to independence

If you are able to get hold of your personal things, you feel better, it ramps up your sense of independence, you seem more in control: and this is a feeling that can accompany residents and patients throughout their daily routines. All kinds of storage units and also easily accessible drawers and compartments help enormously here. Various models with overbed table tops offer additional surfaces and lend themselves as a handy space for reading or dining. The practical and easy-to-use overbed table rounds off this range of furniture perfectly.

+ Advantages

+ Helpful

Storage space, shelves and surfaces that the resident or patient can easily and safely reach.

+ Matching

The décor and functioning of these units match our care beds in every detail.





FMN 204



FMN 205



FMN 206



FMN 210



FMN 211



FMN 212



FMN 215



FMN 250



FMN 220



FMN 222



FMN 230




FMN 231



FMN 232



FMN 233



Optional:

Bedside table with overbed table, rounded; Aluminium tilt fitting at front



Décor selection

Current standard décor selection and trends can be found at fmb-care.de

“In the right kind of surroundings
you can feel better so much more easily.
Together with our customers we design rooms
that create more quality of life.”

GELI FERBER

care & interior with flair

Their rooms in a care home are quite simply the focus of senior citizens' life. And that is certainly true for those in need of a great deal of care. They spend most of their time in their rooms. So for that reason the furnishings should be especially cosy and homely for them. Whereas for the carers, it is vital that equipment and fittings are practical and functional. The care bed and the furniture of the care home should, however, fulfil both kinds of needs.



FINDING THE RIGHT IDEAS TOGETHER

The combination is our strength

These are the keys to good results: teamwork, transfer of knowledge and delight in innovation. With FMB care at your side, you will have the most modern care beds together with the matching room concept and all this from under the one roof.

Know-how and passion


Every kind of property or institution is different, particularly when it comes to interior design. That is why we find it so important that we really get to know the needs and wishes of our customers. Because every establishment, irrespective of whether we are talking

about a single room or a complete building, has its distinctive features, which have to be considered in the planning and implementation of a project. That is why our interior designers work hand in hand with all parties concerned, visualising, if so required, the appropriate, made-to-measure room concepts. From the moment of the first contact right up to the final setting up of the furniture and furnishings at the relevant location, we are there, participating fully with our knowledge and our passion for care products.


+ Advantages

We support the transforming of a resident's room and that of the care establishment by means of diverse planning measures:

A Floor plan



B 3D viewing





A WIDE SELECTION OF DESIGNS

A pleasant ambience




Making a decision for a person to go into care is an exceptionally emotional experience for family members. And many factors certainly play a part in the taking of such a decision. But one crucial issue is the first impression experienced by the family when they first visit the care home. Impressions gained from the moment they go through the entrance until they arrive at the residents' rooms.

Feeling totally at ease

Their rooms in a care home are quite simply the focus of senior citizens' life. And that is certainly true for those in need of a great deal of care. They often spend most of their time in these rooms.

So for that reason the furnishings should be especially cosy and homely for them. For the carers, it is also vital that equipment and fittings are practical and functional. The care bed, the furniture and all the other kinds of furnishings should, however, fulfil both kinds of needs. So when we are developing our Formidabel care beds and also during the design of matching furniture and furnishings, we always incorporate into our work the findings of current studies and, of course, up-to-the-minute trends. Modern wood looks, appealing shapes and attractive fabrics bestow upon a room more cosiness and a discernible feel-good atmosphere.

+ Advantages

-  A harmonious room concept undoubtedly promotes the wellbeing of residents and carers alike.
-  Usefulness and care-friendly features of the furnishings are visually in tune with one another.
-  Family members see this and immediately appreciate the benefits of a cosy environment with homely features.

“As a family-run business we are well aware of the special kind of responsibility needed to understand the requirements of good care and we also realise the responsibility resting on our shoulders to support you with the best solutions in the long term.”

MARTIN FERBER

Together we can **achieve more**

The demands placed on the care sector are changing constantly. Therefore long-term, reliable partnerships are becoming all the more important today. In addition to a comprehensive grasp of the sector, other factors are vital in the success equation: constant dialogue between partners, service offers that provide assistance quickly and directly and also the aspiration of continuously wanting to improve products and services. You will find all these qualities at FMB care and you can rely on them.



FROM EVERY POINT OF VIEW

Understanding

Because of the values that we nurture as a family-run company, here at FMB care we attach huge importance to being your long-term, reliable partner who stands at your side on an equal footing. This philosophy guides our corporate culture, which advocates that at the beginning of a partnership the first thing we do is to let you talk and we listen. Together we then find out what the special requirements are in your project and what you wish to have from us. It makes no difference to us whether you are just interested in accessories, replacements or you need furniture for complete room makeovers. It is always very close to our heart that we are able to support your work in the care sector with the best possible offers and solutions.

Service

Furthermore, when it comes to things like checks, maintenance, servicing and upkeep, FMB care will always be there at your side as your reliable partner. Our care beds are designed and produced for permanent use and to last for many years. Our service offers are on hand to support you and your team so that during this lengthy period virtually no problems will arise. We take care of all manner of coordination tasks and carry out work such as testing and inspections. And we are always available telephonically with our Support Team and our Spare Part Management department to assist you wherever we can. No complications, plenty of understanding.

+ Your direct contact

1 E vertrieb@fmb-care.de
T +49 52 58 / 98 09 018

2 E service@fmb-care.de
T +49 52 58 / 98 09 021



BERLIN CERT

When it comes to the quality of our care beds, we tolerate absolutely no compromises: so as to ensure safety for the users of our beds, we observe all current regulations stipulated for management systems and products found in the field of medical technology.

Independent testing and certification laboratories document our adherence to all requirements and verify our high quality standards.



ISO 13485:2016

This standard attests to comprehensive quality management, representing the design and manufacturing of medical products. The certification obtained in the summer of 2021 states that our customers can rest assured that our care furniture and also our product innovations meet the highest requirements and that our good quality remains constant.



ÖKOPROFIT

FMB care produces care beds completely in tune with the guidelines of Ökoprofit, a German ecological and financial undertaking: aiming at the reduction of operating costs while at the same time conserving natural resources. Our company has been awarded distinctions by the government ministry for our collaboration on sustainable projects focusing on photovoltaic systems, charging stations and on the conversion to LED lighting. FMB care continues to actively commit to environmental protection throughout the company.

SHAPING THE FUTURE

Responsibility

For us as a family-run company, the word responsibility is not just a platitude. On the contrary, it is a part of our corporate DNA. We stand behind all our products and services on offer to you with the greatest degree of conviction.

The person counts

But our perception of responsibility does not stop with our core business. FMB care would like to make an active contribution to improving the situation in the care sector on a permanent basis in order to secure long-term results. Whatever we are developing and whatever the new idea is, we always ask the same question first: How can this benefit people? Every innovation has to make the lives of the residents and patients or those of the carers easier. If we are convinced of this, then we invest with our heart and soul.

Sustainability drives us on

With our products, our production processes and in our business activities, we always remain vigilant about the matter of sustainability. We optimise processes in order to save energy, we use materials whose origin and production conditions are clearly acceptable and we are open to innovations. Moreover, the concept of sustainability plays an important role when it comes to our customer relationships. We are utterly convinced that sustainability and a sense of responsibility are inextricably linked. And that is why we promise to link up with our partners in a long-lasting and sustainable way. A promise that embraces our wish to shape the future together in the best possible way.


+ How to reach us

FMB care
Am Schlinge 8
D-33154 Salzkotten

E info@fmb-care.de
T +49 52 58 / 98 09 00

www.fmb-care.de

Notes





FMB care GmbH

Am Schlinge 8

D-33154 Salzkotten

+49 (0) 52 58 / 98 09 00

info@fmb-care.de

fmb-care.de